

## Relative Activity of OPTIZYME™ Restriction Enzymes in 1X Buffers

The 10X Buffer supplied with each restriction enzyme is optimized to give 100% activity at the 1X concentration. This table may be used to select the best buffer for digestions with multiple restriction enzymes. Enzyme activity is expressed as a percent of the activity obtained with the optimized buffer for each enzyme in a one-hour digest.

Part Number	OPTIZYME* Restriction Enzyme	Recommended buffer for 100% activity	Enzyme activity in 1X OPTIZYME* buffers, %					Optimal reaction temperature	Activity at 37° C, %
			Buffer 1	Buffer 2	Buffer 3	Buffer 4	Buffer 5		
BP8069	AarI	Unique	NR (+oligo)	NR (+oligo)	0-20 (+oligo)	NR (+oligo)	0-20 (+oligo)	37°C	100%
BP8041	AatII	Buffer 4	50-100	20-50	0-20	100	0-20	37°C	100%
BP8075	Alol	Buffer 5	0-20	0-20	0-20	20-50	100	30°C	20%
BP8015	Alul	Buffer 4	50-100	0-20	0-20	100	0-20	37°C	100%
BP8059	Alw44I (ApaLI)	Buffer 4	50-100	100	0-20	100	50-100	37°C	100%
BP8025	Apal	Buffer 1	100	20-50	0-20	20-50	0-20	37°C	100%
BP8035	AvaI (Eco88I)	Buffer 4	100	50-100	0-20	100	0-20	37°C	100%
BP8043	AvaII (Eco47I)	Buffer 5	0-20	50-100	50-100	50-100	100	37°C	100%
BP8039	BalI (MscI)	Buffer 5	0-20	20-50	0-20	20-50	100	37°C	100%
BP8005	BamHI	Unique	20-50*	100	20-50	100*	50-100*	37°C	100%
BP8053	BclI	Buffer 2	20-50	100	20-50	100*	20-50	55°C	50%
BP8046	BglI	Buffer 3	0-20	50-100	100	0-20	100	37°C	100%
BP8014	BglII	Buffer 3	0-20	20-50	100	0-20	50-100	37°C	100%
BP8072	BpiI (BbsI)	Buffer 2	20-50	100	50-100	50-100	50-100	37°C	100%
BP8071	Bsh1236I (BstUI)	Buffer 5	0-20	0-20	50-100	20-50	100	37°C	100%
BP8078	BshTI (AgeI)	Buffer 3	0-20	20-50	100	20-50	50-100	37°C	100%
BP8036	BssHII (PauI)	Buffer 5	0-20	0-20	100	0-20	100	37°C	100%
BP8038	BstEII (Eco91I)	Buffer 3	20-50	20-50	100	NR	50-100	37°C	100%
BP8081	Cfr9I (XmaI)	Unique	0-20	0-20	0-20	20-50	0-20	37°C	100%
BP8024	Clal	Buffer 4	20-50	20-50	20-50	100	20-50	37°C	100%
BP8068	Csp6I (CviQI)	Buffer 1	100	50-100	0-20	50-100	0-20	37°C	100%
BP8060	DdeI (HpyF31)	Buffer 4	20-50	20-50	20-50	100	20-50	37°C	100%
BP8009	DpnI	Buffer 4	100	100	50-100	100	50-100	37°C	100%
BP8026	DraI	Buffer 4	50-100	50-100	20-50	100	20-50	37°C	100%
BP8080	Ecl136II	Unique	50-100	20-50	0-20	50-100	0-20	37°C	100%
BP8066	Eco57I (AclI)	Buffer 2 +SAM	100 (+SAM)	100 (+SAM)	20-50 (+SAM)	50-100 (+SAM)	20-50 (+SAM)	37°C	100%

## Relative Activity of OPTIZYME™ Restriction Enzymes in 1X Buffers

The 10X Buffer supplied with each restriction enzyme is optimized to give 100% activity at the 1X concentration. This table may be used to select the best buffer for digestions with multiple restriction enzymes. Enzyme activity is expressed as a percent of the activity obtained with the optimized buffer for each enzyme in a one-hour digest.

Part Number	OPTIZYME* Restriction Enzyme	Recommended buffer for 100% activity	Enzyme activity in 1X OPTIZYME* buffers, %					Optimal reaction temperature	Activity at 37° C, %
			Buffer 1	Buffer 2	Buffer 3	Buffer 4	Buffer 5		
BP8003	EcoRI	Unique	0-20	NR	100	NR	100*	37°C	100%
BP8054	EcoRI (HC)	Unique	0-20	NR	100	NR	100*	37°C	100%
BP8012	EcoRV (Eco32I)	Buffer 5	0-20	50-100	50-100	20-50	100	37°C	100%
BP8070	Esp3I (BsmBI)	Buffer 4	100 (+DTT)	20-50 (+DTT)	0-20 (+DTT)	100 (+DDT)	0-20 (+DTT)	37°C	100%
BP8002	HaeIII (BsuRI)	Buffer 5	0-20	0-20	0-20	50-100	100	37°C	100%
BP8034	HincII (HindII)	Buffer 4	50-100	50-100	20-50	100	50-100	37°C	100%
BP8006	HindIII	Buffer 5	0-20	20-50	0-20	50-100	100	37°C	100%
BP8051	Hinfl	Buffer 5	0-20	20-50	50-100	50-100	100	37°C	100%
BP8049	HpaI (KspAI)	Buffer 1	100	50-100*	20-50	100*	20-50	37°C	100%
BP8032	HpaII (MspI)	Buffer 4	50-100	50-100	0-20	100	20-50	37°C	100%
BP8079	Hpy8I (MjaIV)	Buffer 4	50-100	50-100	0-20	100	20-50	37°C	100%
BP8083	KpnI	Unique	20-50	0-20	0-20	20-50	0-20	37°C	100%
BP8067	LguI (SapI)	Buffer 4	20-50	50-100	20-50	100	20-50	37°C	100%
BP8021	MluI	Buffer 5	0-20	20-50	50-100	20-50	100	37°C	100%
BP8048	MspI (HpaII)	Buffer 4	50-100	50-100	0-20	100	0-20	37°C	100%
BP8057	NaeI	Buffer 4	50-100	20-50	0-20	100	0-20	37°C	100%
BP8017	NcoI	Buffer 4	20-50	20-50	20-50	100	50-100	37°C	100%
BP8020	NdeI	Buffer 3	0-20	0-20	100	0-20	50-100	37°C	100%
BP8019	NheI	Buffer 4	100	20-50	0-20	100	0-20	37°C	100%
BP8004	NotI	Buffer 3	0-20	20-50	100	0-20	20-50	37°C	100%
BP8058	NsiI (Mph1103I)	Buffer 5	0-20	50-100	20-50	50-100	100	37°C	100%
BP8073	PasI	Unique	NR	NR	NR	NR	NR	55°C	30%
BP8077	PfoI	Buffer 4	0-20	20-50	50-100	100	0-20	37°C	100%
BP8001	PstI	Buffer 3	50-100	50-100	100	50-100	100	37°C	100%
BP8050	PvuI	Buffer 5	0-20	20-50	50-100	50-100	100	37°C	100%
BP8022	PvuII	Buffer 2	50-100*	100	20-50	20-50*	50-100	37°C	100%

## Relative Activity of OPTIZYME™ Restriction Enzymes in 1X Buffers

The 10X Buffer supplied with each restriction enzyme is optimized to give 100% activity at the 1X concentration. This table may be used to select the best buffer for digestions with multiple restriction enzymes. Enzyme activity is expressed as a percent of the activity obtained with the optimized buffer for each enzyme in a one-hour digest.

Part Number	OPTIZYME* Restriction Enzyme	Recommended buffer for 100% activity	Enzyme activity in 1X OPTIZYME* buffers, %					Optimal reaction temperature	Activity at 37° C, %
			Buffer 1	Buffer 2	Buffer 3	Buffer 4	Buffer 5		
BP8000	RsaI	Buffer 4	50-100	20-50	0-20	100	0-20	37°C	100%
BP8016	SacI	Unique	50-100	20-50	0-20	50-100	0-20	37°C	100%
BP8023	SacII (Cfr42I)	Buffer 1	100	50-100	0-20	50-100	0-20	37°C	100%
BP8013	Sall	Buffer 3	0-20	0-20	100	0-20	20-50	37°C	100%
BP8030	Sau3AI (Bsp143I)	Unique	20-50	20-50	0-20	50-100	0-20	37°C	100%
BP8037	Scal	Unique	0-20	0-20	0-20	0-20	0-20	37°C	100%
BP8076	SfaAI (AsiSI)	Buffer 4	50-100	0-20	0-20	100	0-20	37°C	100%
BP8011	SmaI	Buffer 4	50-100	0-20	0-20	100	0-20	30°C	50%
BP8018	SpeI (BcuI)	Buffer 4	50-100	50-100	0-20	100	20-50	37°C	100%
BP8029	SphI (PaeI)	Buffer 1	100	50-100	0-20	50-100	0-20	37°C	100%
BP8027	StuI (Eco147I)	Buffer 1	100	50-100	20-50	50-100	20-50	37°C	100%
BP8007	TaqI	Unique	0-20	20-50	20-50	20-50	20-50	65°C	10%
BP8064	Tru9I (MseI)	Buffer 5	50-100	50-100	20-50	50-100	100	65°C	10%
BP8055	VspI (AseI)	Buffer 3	0-20	50-100	100	100	20-50	37°C	100%
BP8008	XbaI	Buffer 4	50-100	50-100	20-50	100	0-20	37°C	100%
BP8010	XhoI	Buffer 5	0-20	50-100	50-100	20-50	100	37°C	100%
BP8082	XmaII (AvrII)	Buffer 4	20-50	50-100	50-100	100	50-100	37°C	100%
BP8052	XmnI (PdmI)	Buffer 4	20-50	50-100	0-20	100	0-20	37°C	100%

### Note

\* Star activity appears at a greater than 5-fold overdigestion (5 units x 1 hour).

NR – buffer is not recommended, because of high star activity.