

Save up to 25% on your manual western blot workflow

Start the year with great savings on Thermo Scientific™ manual western blot products. Choose from membranes, blocking buffers, western blot signal enhancers and secondary antibodies.

To claim your discount, simply place your order with Fisher Scientific. This offer expires on 30 June 2016. Please quote promotion codes when ordering:

25% off - Quote RCS011

15% off - Quote RCS012

Membranes for western blotting

To start your western blot, choose from nitrocellulose and PVDF transfer membranes in pre-cut sheets.

25% off

Description	Catalogue No.	Size
PVDF Transfer Membrane, 0.45 µm, 10 cm x 10 cm	11804131	10 sheets
Nitrocellulose Membrane, 0.2 µm, 8 cm x 8 cm	10096363	15 sheets

Blocking buffers for western blotting

Block unreacted sites on the membrane to reduce the amount of nonspecific binding. To achieve the highest signal-to-noise ratio possible for your blots, choose between two options: Thermo Scientific™ StartingBlock™ Buffers can be used for ELISA, western blot, dot blot and immunohisto-chemistry. Thermo Scientific™ SEA BLOCK Blocking Buffer can be used for ELISA, western blot and dot blot.

25% off

Description	Catalogue no.	Size
StartingBlock (PBS) Blocking Buffer	10108313	1 L
StartingBlock (TBS) Blocking Buffer	10371804	1 L
SEA BLOCK Blocking Buffer	10782925	500 mL

Western blot signal enhancers

When a protein or antigen is difficult to detect because of low abundance or poor immunoreactivity, use of Thermo Scientific™ SuperSignal™ Western Blot Enhancer can significantly reduce background and enhance detection of low-abundance and weakly immunoreactive antigens.

25% off

Description	Catalogue no.	Size
SuperSignal Western Blot Enhancer	10157958	500 mL

Filter paper for western blotting

Sheets of blotting paper are necessary components of transfer sandwiches and cassettes that typically must be assembled for various kinds of protein or nucleic acid, gel-to-membrane, transfer protocols.

Thermo Scientific™ Western Blotting Filter Papers are pre-cut cotton sheets for wet or semi-dry, passive or electrophoretic transfer of proteins from polyacrylamide gels (SDS-PAGE) to PVDF, nitrocellulose, or other membranes.

15% off

Description	Catalogue no.	Size
Western Blotting Filter Paper, 8 cm x 13.5 cm	15115214	100 sheets
Western Blotting Filter Paper, 7 cm x 8.4 cm	15105214	100 sheets

Thermo Scientific™ Power Blotter

For high-performance semi-dry blotting, use the Thermo Scientific™ Power Blotter. The unit is designed to accomplish semi-dry transfer of proteins (10–300 kDa) from polyacrylamide gels to nitrocellulose or PVDF membranes in just 5–10 minutes. The integrated power supply, easy-touch programming and versatility to accommodate up to four mini gels or two midi gels at a time, will ensure you achieve rapid and efficient protein transfer every time.

25% off

Description	Catalogue no.	Size
Power Blotter Welcome Pack – UK Plug	15309134	1 unit
Power Blotter Welcome Pack – EU Plug	15319134	1 unit

Secondary antibodies for western blotting

You need the right secondary antibodies to get clean and clear results. Our secondary antibodies, conjugated with HRP or Infra-red dyes, provide the specificity and sensitivity to help you get the Western-blot data you need to advance your research. They are also fully validated, eliminating the need to screen numerous antibodies to find the correct one.

25% off

Description	Catalogue no.	Size
Goat anti-Mouse IgG (H+L) Secondary Antibody, HRP conjugate	10158113	2 mL
Goat anti-Rabbit IgG (H+L) Secondary Antibody, HRP conjugate	11859140	2 mL
Goat anti-Mouse IgG (H+L) Cross Adsorbed Secondary Antibody, DyLight 680 conjugate	10797775	500 µL
Goat anti-Mouse IgG (H+L) Cross Adsorbed Secondary Antibody, DyLight 800 conjugate	13477327	500 µg
Rabbit anti-Mouse IgG (H+L) Cross Adsorbed Secondary Antibody, DyLight 680 conjugate	13437317	500 µg
Rabbit anti-Mouse IgG (H+L) Cross Adsorbed Secondary Antibody, DyLight 800 conjugate	13457317	500 µg

For more information and to order, please contact Fisher Scientific



ThermoFisher
SCIENTIFIC

Austria: +43 (0)800 20 88 40 **Belgium:** +32 (0)056 260 260 **Denmark:** +45 70 27 99 20
Germany: +49 2304 932-5 **Ireland:** +353 (0)1 885 5854 **Italy:** +39 02 950 59478
Finland: +358 (0)9 8027 6280 **France:** +33 (0)3 88 67 14 14 **Netherlands:** +31 (0)20 487 70 00
Norway: +47 22 95 59 59 **Portugal:** +351 21 425 33 50 **Spain:** +34 902 239 303
Sweden: +46 31 352 32 00 **Switzerland:** +41 (0)56 618 41 11 **UK:** +44 (0)1509 555 500

* Terms and Conditions: These promotions are open to customers in Austria, Belgium, Denmark, Finland, France, Germany, Ireland, Italy, Netherlands, Norway, Portugal, Spain, Sweden, Switzerland and UK. To be eligible for the promotion price, customers must quote the relevant promotion code when ordering. Promotion is valid until 30 June 2016 or until promotional supplies are depleted, whichever comes first. Cannot be combined with any other promotions or discounts. Void where prohibited by federal, state, provincial, or local laws or regulation or agency/institutional policy.

For Research Use Only. Not for use in diagnostic procedures. © 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. CO37923 0116