

All laboratory labels for BMP71



Brady offers durable identification labels for any laboratory sample. The labels are engineered to stay attached and remain legible in the most challenging laboratory conditions when printed with the BMP71 Label Printer.

Cryo Vial Identification



B-490, B-492

FreezerBondz™ ultra-thin polyester label to label both frozen and room temperature vials, tubes and canes. Superior chemical and abrasion resistance. B-490 is designed for flat surfaces, B-492 for round surfaces.

Withstands







Well Plate Identification



B-430

Clear polyester label with excellent chemical resistance. Ideal for well plate identification. Applies well to a wide range of surfaces.

Withstands





Tube, Vial & Straw Identification



Ultra-thin polyester label with clear tail that protects print and allows liquid levels to be viewed.

Withstands











B-427

Self-laminating vinyl label with excellent abrasion and smudge resistance suitable for small surface areas of straws and rods, withstanding long-term storage.

Withstands









Bottle tamper seal labels



B-351

Vinyl label designed to fracture easily in order to show signs of sample tampering and to prevent one-piece label removal. Provides tracking information with the added feature of security.

Withstands



Freezer

Dish/Bottle Identification



Thick nylon cloth label that is easy to apply and conforms well to various shapes for excellent adhesion.

Withstands







Tissue Cassette Identification



B-498

Vinyl cloth label to apply directly on tissue cassettes. Can also be easily applied on pure paraffin and on top of older labels.

Withstands





B-481

Stainerbondz polyester label with extreme chemical resistance. Performs great in slide staining processes.

Withstands



Н

B-488

Polyester label with excellent chemical resistance. Ideal to identify slides and other flat surfaces.

Withstands





General Purpose Identification



TONGUE DEPRESSORS

B-423

Highly adhesive white label for drawers. shelves, cabinets and general labelling.

Withstands







B-595

Highly adhesive label for drawers, shelves, cabinets and general labelling. Available in white, red, yellow, green, blue, orange, black & clear.

Withstands



Freeze

B-424

Low cost, general purpose labeling.

Withstands

