



SuperClear[®] Specimen Collection and Transport Tubes

Complete Safety in Specimen Shipment

✓ 95kPa certified as designated by the IATA and U.S. DOT Standard

Find it at fishersci.eu

 **fisher scientific**
part of Thermo Fisher Scientific

SuperClear® Specimen Collection and Transport Tubes


Complete Safety In Specimen Shipment

Labcon's SuperClear® Specimen Collection and Transport Tubes provide certified safety when storing and shipping biological material. These leakproof freestanding tubes meet the IATA and U.S. DOT 95kPa performance specifications legally required when shipping samples by ground or air. This ensures transport tubes are leakproof under differential pressure conditions and extreme temperature shifts as outlined in IATA DGR UN3373 and U.S. DOT HMR, 49 C.F.R., Parts 171-180.

Made from 100% virgin USP Class-VI medical-grade resins, SuperClear® Specimen Collection and Transport Tubes are suitable for use in medical, diagnostic and research labs and are certified contamination-free as mandated by the U.S. CDC. These certified sterile tubes have a Sterility Assurance Level (SAL) of 10^{-6} , as required for diagnostic assays. Provides certified safety when storing and shipping biological material.

- ✓ **95kPa certified**
 Meets legal requirements for IATA/U.S.DOT/UN3373 air & ground shipping
- ✓ **Leakproof**
 Certified leakproof under differential pressure & extreme temperature shifts
- ✓ **Validated Sterility**
 10^{-6} SAL, as required for diagnostic & medical assays

Product Features



Quality & Compliance



Packaging	Sterile	Qty Pk	Qty Cs	Cat. No.
5 mL SuperClear® Specimen Collection and Transport Tubes				
Individually Wrapped	✓	1	100	16696846
Bag		50	500	16606856
Bag	✓	50	500	16616856
Bulk		500	1,000	16626856

Packaging	Sterile	Qty Pk	Qty Cs	Cat. No.
10 mL SuperClear® Specimen Collection and Transport Tubes				
Individually Wrapped	✓	1	100	16636856
Bag		50	500	16646856
Bag	✓	50	500	16656856
Bulk		500	1,000	16666856

Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

00-0000 BG 09/20 BNXXXXXXXX

