

Thermo Scientific Direct Heat CO2 Incubators Comparison Chart

Product	Capacity (cu ft)	In Chamber Class 100 HEPA Air Filtration	High Heat Sterilization Cycle	External Humidity Source & Control	100% Solid Copper Interior
Heracell 150i and 240i with Copper Chambers Interactive touch screen simplicity with 100% solid copper chamber	5.4 or 8.4		X		X
Forma* Steri-Cult* CO2 Incubators Ultimate Protection, control, and capacity	8.2 or 11.4	X	X	X	Optional shelves and pan
Forma* 310 Direct-Heat CO2 Incubators Superior value, basic culturing	6.5	Optional			Optional shelves and pan
Forma* Steri-Cycle* CO2 Incubators Patented HEPA filtration and sterilization	6.5	X	X		Optional shelves and pan
Heracell* 150i and 240i CO2 Incubators with Stainless-Steel Chambers Interactive touch screen simplicity	5.4 or 8.4		X		